



Be Right™



Flask, Erlenmeyer, polymethylpentene capacity 250 mL

Kat. číslo produktu:	2089846
EUR Cena bez DPH:	Kontaktujte nás
Dostupné	

Polymethylpentene flask with polypropylene screwcap features clarity, chemical resistance, strength, and autoclavability (Do not engage threads during cycle). Useful for titrations, preparation and storage of media.

Technické údaje

Jednotky:	1 each
Kapacita:	250 mL
Materiál:	Polymethylpentene
Typ:	Erlenmeyer