



Glass Fiber Filter, Diameter 70 mm, Preweighed for Standard Methods TSS procedure, 100/pk

Kat. číslo produktu: 2546135
EUR Cena bez DPH: Kontaktujte nás
Dátum odoslania neuvedené

Individually preweighed filter papers, preweighed in compliance with Standard Methods Total Suspended Solids procedure. Each in a planchet with tare weight listed. Each tare weight meets Standard Methods specifications.

Technické údaje

Jednotky:	100/pk
Materiál:	Glass Microfiber
Priemer:	70 mm
Veľkosť balenia:	100 /pk
Veľkosť póru:	1.5 µm