



Test kit, phosphorus, orthophosphate (reactive) , model PO-19

Kat. číslo produktu:

224800

Nebezpečné

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené



Method / Chemistry: Color disc / Ascorbic acid. Range (mg/L): 0-1, Smallest Increment- 0.02. Range (mg/L): 0-5, Smallest Increment- 0.1. Range (mg/L): 0-50, Smallest Increment- 1. Approx. Number of Tests: 100. Case Style: 12x7.25x6 inches Blue Polypropylene Carrying Case.

Results in less than 5 minutes for clean water samples

Utilizes simple color disc comparator for accurate measurements

Model PO-19A contains filtration apparatus for turbid samples

Technické údaje

Model:	PO-19
Model testovacej súpravy:	PO-19
Názov metódy:	Color disc/Ascorbic acid
Parameter:	Phosphorus
Platforma:	Farebný disk
Počet testov:	100
Prepravná hmotnosť:	3,1
Rozsah merania:	0-0,8 mg/L PO ₄
Rozsah merania (2):	0-4 mg/L PO ₄
Smallest increment:	0.02
	0.10
	1
Štandardný rozsah merania:	0 - 0.8 mg/L PO ₄ ; 0 - 4 mg/L PO ₄ ; 0 - 40 mg/L PO ₄
Štýl kufríka:	R