



Test kit, phosphorus, orthophosphate (reactive) , model PO-19A

Kat. číslo produktu:

224801

Nebezpečné

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené



Simple, low-cost color disc method is the most accurate visual test kit method

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Technické údaje

Model: PO-19A

Model testovacej súpravy: PO-19A

(Includes filtration for turbid samples.)

Názov metódy: Color disc/Ascorbic acid

Parameter: Phosphorus

Platforma: Farebný disk

Počet testov: 100

Prepravná hmotnosť: 3,7

Rozsah merania: 0-0,8 mg/L PO₄

Rozsah merania (2): 0-4 mg/L PO₄

Smallest increment: 0.02

Štandardný rozsah merania: 0 - 0.8 mg/L PO₄; 0 - 4 mg/L PO₄; 0 - 40 mg/L PO₄

Štýl kufrika: F