



## Test kit, phosphate, total ortho-/meta-, model PO-24

**Kat. číslo produktu:**

**225001**

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené

Nebezpečné



Kit provides measurement of ortho- (reactive) and meta- (condensed) phosphate (requires hot acid digestion).

**Kit provides rapid measurement of ortho- (reactive) and meta- (condensed) phosphate (requires hot acid digestion).**

---

### Technické údaje

|                            |                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------|
| Model:                     | PO-24                                                                                   |
| Model testovacej súpravy:  | PO-24                                                                                   |
| Názov metódy:              | Color disc/Ascorbic acid                                                                |
| Parameter:                 | Phosphate                                                                               |
| Platforma:                 | Farebný disk                                                                            |
| Počet testov:              | 100 čirá                                                                                |
|                            | 100 zakalená                                                                            |
| Prepravná hmotnosť:        | 4,5                                                                                     |
| Rozsah merania:            | 0-0,8 PO <sub>4</sub>                                                                   |
| Rozsah merania (2):        | 0-4 mg/L PO <sub>4</sub>                                                                |
| Smallest increment:        | 0.02                                                                                    |
|                            | 0.1                                                                                     |
|                            | 1                                                                                       |
| Štandardný rozsah merania: | 0 - 0.8 mg/L PO <sub>4</sub> ; 0 - 4 mg/L PO <sub>4</sub> ; 0 - 40 mg/L PO <sub>4</sub> |
| Štýl kufrika:              | F                                                                                       |