



Be Right™



## Testkit, glycol visual, model EG-1

Kat. číslo produktu:

2186400

Nebezpečné

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené



Semiquantitative determination of glycol in water or oil. Visual appearance of colour indicates glycol presence.

Ethylene glycol (anti-freeze) contamination of automotive lubricants is an indication of a cracked engine block. Glycol presence also can indicate leaks in cooling systems. Test Kit Model EG-1 presents a purple color if ethylene glycol and other 1,2 glycols are present in oil or water. The kit contains approximately 25 tests and comes with a 14.4" x 8" x 7" Blue Polypropylene Carrying Case.

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Technické údaje

Model testovacej súpravy:	EG-1
Názov metódy:	Visual/ Appearance of color
Parameter:	Glycol - Ethylene glycol in oil or water
Počet testov:	25
Prepravná hmotnosť:	3,5
Rozsah merania:	Přítomnost nebo nepřítomnost glykolu, 0- více než 150 mg/L.
Štýl kufrička:	G