



Total Hardness Test Kit, Model 5-EP

Kat. číslo produktu:

145400

Nebezpečné

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené



Water hardness is caused almost entirely by calcium and magnesium ion. Other di- and trivalent metals have a similar effect, but usually are not present in high enough concentration in potable waters to cause problems. Hardness increases soap consumption in laundries and causes scale in boilers. With the Model 5-EP test kit, simply count the number of drops required to change the solution color to determine hardness concentration. Utilizes drop count titration with a range of 1-30 gpg. Contains approximately 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Model testovacej súpravy:	5-EP
Názov metódy:	Drop count titration/ EDTA-powder pillows
Parameter:	Hardness, Total - As CaCO_3
Platforma:	Titrácia s počítaním kvapiek
Počet testov:	100
Rozsah merania:	1 - 30 gpg CaCO_3
Smallest increment:	1 gpg
Štýl kufríka:	B