



Manganese test kit, model MN-5

Kat. číslo produktu:

146700

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené

Nebezpečné



Manganese causes a bitter taste in water, and at concentrations above 0.1 mg/L, it causes objectionable stains on laundry and plumbing fixtures. Manganese control also is important in the beverage, paper, and textile industries, and in dye production and food processing. An economical solution, this kit uses the simple, accurate color disc method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Technické údaje

Model:	MN-5
Model testovacej súpravy:	MN-5
Názov metódy:	Color disc/Cold periodate
Parameter:	Manganese
Platforma:	Farebný disk
Počet testov:	100
Rozsah merania:	0 - 3,0 mg/L Mn
Smallest increment:	0.1
Štandardný rozsah merania:	0 - 3.0 mg/L Mn
Štýl kufríka:	D