



Be Right™



## Test kit, hypochlorite, model CN-HRDT

Kat. číslo produkta:

2687100

Nebezpečné

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené



Hypochlorite Kit, Model CN-HRDT with Digital Titrator

Utilizing the Hach Digital Titrator, this kit is used to determine hypochlorite, high range (Bleach) solution concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 5 - 15% as Cl<sub>2</sub>: 50 - 150 g/L as Cl<sub>2</sub>. Contains 100 tests.

**Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration**

**Rugged titration device is ideal for field use, with no glass burets required**

**Designed for on-site analysis of bulk hypochlorite deliveries or bulk storage solutions at drinking water or chloramination facilities**

**Verify concentration of hypochlorite solutions at time of delivery**

**Kit contains all required reagents, solutions, and glassware**

### Technické údaje

Číslo metódy:	10100
Model testovacej súpravy:	CN-HRDT
Názov metódy:	Digital Titrator/
	Thiosulfate
Parameter:	Hypochlorite, high range (Bleach) - As Cl <sub>2</sub>
Platforma:	Digitálny titrátor
Počet testov:	100
Rozsah merania:	5 - 15 % = 50-150 g/L Cl <sub>2</sub>
Smallest increment:	0.05%
Štýl kufríka:	G