



*Be Right*TM



Test kit, hypochlorite, model CN-HR

Kat. číslo produkta:

2687200

Nebezpečné

EUR Cena bez DPH

Kontaktujte nás

Dátum odoslania neuvedené

三



Hypochlorite Kit, Model CN-HR - Drop Count

Utilizing the simple, economical drop count titration method, this kit is used to determine high range hypochlorite solution concentrations (bleach solutions). Range: 5-15%, contains 100 tests.

Simple, economical titration method

Kit contains all required reagents, solutions, and glassware

Technické údaje

Model testovacej súpravy:

CN-HR

Názov metódy:

Drop count titration/

Thiosulfate

Hypochlorite, high range - (Bleach) As Cl.

Parameter:

Drop Count Titration

Počet testov:

100

Bozsah moranja:

50 - 150 g/L Cl_x (5-15%)

Smallest increment:

0.5%