



Cylinder, graduated, 25 mL $\pm$ 0.3 mL, 0.5 mL divisions (white markings)

Kat. číslo produktu: 50840
EUR Cena bez DPH: Kontaktujte nás
Dátum odoslania neuvedené

Calibrated "to deliver" glass cylinder. Single metric scale in durable white markings. Bumper guards included for all sizes $\geq$ 25 mL.

Technické údaje

Certifikácia:	Class B
Kapacita:	25 mL
Materiál skrine:	Glass
Metóda merania:	To Contain
Typ:	S jednou stupnicou
Veľkosť balenia:	1 each